## ExamVision Magnifications – Galilean System



Please note that the actual perceived magnification can vary, depending on working distance and prescription.

These examples are to be used as a guide and represent angular magnification, based on a working distance of 45 cm.



No magnification

## Galilean HD



Magnification	2.3x	2.8x	3.3x	3.8x
Field of View	12 cm	10.5 cm	8 cm	6.5 cm
Depth of Field up to	28 cm	20 cm	15 cm	10 cm

## **Galilean Essential**



Ma	agnification	2.5x	3.0x
Fie	eld of View	11 cm	8.3 cm
De	epth of Field up to	25 cm	15 cm